Rev 3.0

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- 1 Centers for Disease Control and Prevention. Interim Clinical Guidance for Management of Patients with Confirmed Coronavirus Disease (COVID-19). Retrieved from <https://www.cdc.gov/coronavirus/2019-ncov/hcp/clinical-guidance-management-patients.html> March 22, 2020
- 2 WHO. Clinical management of severe acute respiratory infection (SARI) when COVID-19 disease is suspected: interim guidance, 13 March 2020. Retrieved from <https://apps.who.int/iris/bitstream/handle/10665/331446/WHO-2019-nCoV-clinical-2020.4-eng.pdf?sequence=1&isAllowed=y> March 17, 2020.
- 3 Guan WJ, Ni ZY, Hu Y, et al. Clinical Characteristics of Coronavirus Disease 2019 in China [published online ahead of print, 2020 Feb 28]. *N Engl J Med.* 2020
- 4 Wang D, Hu B, Hu C, et al. Clinical Characteristics of 138 Hospitalized Patients With 2019 Novel Coronavirus-Infected Pneumonia in Wuhan, China [published online ahead of print, 2020 Feb 7]. *JAMA.* 2020
- 5 Murthy S, Gomersall CD, Fowler RA. Care for Critically Ill Patients With COVID-19 [published online ahead of print, 2020 Mar 11]. *JAMA.* 2020
- 6 NHS. Clinical guide for the use of acute non-invasive ventilation in adult patients hospitalised with suspected or confirmed coronavirus during the coronavirus pandemic [19 March 2020, Version 1]. Retrieved from <https://www.england.nhs.uk/coronavirus/wpcontent/uploads/sites/52/2020/03/clinical-guide-acute-niv-ventilation-v1-19-march-2020.pdf> March 22, 2020.
- 7 Handbook of COVID-19 Prevention and Treatment. The First Affiliated Hospital, Zhejiang University School of Medicine Compiled According to Clinical Experience. Retrieved from <https://covid-19.alibabacloud.com/> on March 18, 2020.
- 8 Xie, J., Tong, Z., Guan, X. et al. Critical care crisis and some recommendations during the COVID-19 epidemic in China. *Intensive Care Med* (2020). <https://doi.org/10.1007/s00134-020-05979-7>
- 9 Hick, JL, Hanfling D, Wynia MK, et al. Duty to Plan: Health Care, Crisis Standards of Care, and Novel Coronavirus SARS-CoV-2. *NAM Perspectives.* Discussion paper. National Academy of Medicine. Washington, DC. 2020. Retrieved from <https://doi.org/10.31478/202003b> March 17, 2020.
- 10 NHS. Guidance for the role and use of non-invasive respiratory support in adult patients with coronavirus (confirmed or suspected) [26 March 2020 Version 2].
- 11 U.S. Food and Drug Administration (FDA). Ventilator Supply Mitigation Strategies: Letter to Health Care Providers. Retrieved from <https://www.fda.gov/medical-devices/letters-health-care-providers/ventilator-supply-mitigation-strategies-letter-health-care-providers> March 30, 2020.
- 12 Italian Thoracic Society (ITS-AITO), Association for the Rehabilitation of Respiratory Failure (ARIR). Joint statement on the role of respiratory rehabilitation in the COVID-19 crisis: the Italian position paper. March 8, 2020.
- 13 Medicines & Healthcare Products Regulatory Agency (MHRA). Rapidly Manufactured Ventilator Systems (RMVS). Retrieved from https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/874279/RMVS001_Rapidly_Manufactured_Ventilator_Specification__PDF.pdf March 23, 2020
- 14 US Department of Health and Human Services. Food and Drug Administration. Center for Devices and Radiological Health. Enforcement Policy for Ventilators and Accessories and Other Respiratory Devices During the Coronavirus Disease 2019 (COVID-19) Public Health Emergency. Guidance for Industry and Food and Drug Administration Staff. March 2020

- 15 Robinson K, Polos P, Stewart JJ. Intubated BiPAP Use to Delay Mechanical Ventilation in COVID-19. Retrieved from <https://www.syneoshealth.com/sites/default/files/documents/Intubated%20BiPAP%20Use%20to%20Delay%20Mechanical%20Ventilation%20in%20COVID-19.pdf> March 30, 2020
- 16 American Association for Respiratory Care. Joint Statement on Multiple Patients Per Ventilator. Retrieved from https://www.aarc.org/wpcontent/uploads/2020/03/032620-COVID-19-press-release.pdf?utm_source=informz&utm_medium=email&utm_campaign=aarcprotech&_zs=M4vKN1&_zl=f3if5 March 28, 2020.
- 17 Tran K, Cimon K, Severn M, Pessoa-Silva CL, Conly J. Aerosol generating procedures and risk of transmission of acute respiratory infections to healthcare workers: a systematic review. *PLoS One*. 2012;7(4):e35797. doi:10.1371/journal.pone.0035797
- 18 Alhazzani W, Møller MH, Arabi YM, et al. Surviving Sepsis Campaign: Guidelines on the Management of Critically Ill Adults with Coronavirus Disease 2019 (COVID-19). *European Society of Intensive Care Medicine and the Society of Critical Care Medicine*. 2020.
- 19 Simonds AK, Hanak A, Chatwin M, et al. Evaluation of droplet dispersion during non-invasive ventilation, oxygen therapy, nebuliser treatment and chest physiotherapy in clinical practice: implications for management of pandemic influenza and other airborne infections. *Health Technol Assess*. 2010
- 20 Cascella M, Rajnik M, Cuomo A, et al. Features, Evaluation and Treatment Coronavirus (COVID-19) [Updated 2020 Mar 8]. In: *StatPearls* [Internet]. Treasure Island (FL): StatPearls Publishing; 2020 Jan-. Available from: <https://www.ncbi.nlm.nih.gov/books/NBK554776/>
- 21 Gonzalez J, Maisonobe J, Oranger M, Mendoza-Ruiz A. Homecare respiratory equipment for patients suspected to be infected with the respiratory virus Covid19. Retrieved from <http://splf.fr/wp-content/uploads/2020/03/LES-PROCEDURES-DU-GAVO2-ProtectionVirale2020-J14mars2020.pdf> March 18, 2020
- 22 CDC. Interim Infection Prevention and Control Recommendations for Patients with Suspected or Confirmed Coronavirus Disease 2019 (COVID-19) in Healthcare Settings. Retrieved from <https://www.cdc.gov/coronavirus/2019-ncov/infection-control/control-recommendations.html> March 17, 2020.
- 23 FDA US Food & Drug Administration. Emergency Use Authorization Letter, March 24, 2020. Retrieved from <https://fda.gov/media/136423/download>.
- 24 FDA US Food & Drug Administration. Authorized Ventilators, Ventilator Tubing Connectors, and Ventilations Accessories. Retrieved from <https://www.fda.gov/media/136528/download>
- 25 Lumis and AirCurve are essentially the same device, with different country variant naming, specific indications for use of variants of the Lumis and AirCurve device may vary by region.
- 26 Environmental Protection Agency. List N: disinfectants for use against SARS-CoV-2. EPA website. Retrieved from <https://www.epa.gov/pesticide-registration/list-n-disinfectants-use-against-sars-cov-2> March 22, 2020.
- 27 Australian Government Department of Health. Environmental cleaning and disinfecting principles – Version 1 (10/03/2020); Coronavirus disease (COVID-19). Retrieved from <https://www.health.gov.au/sites/default/files/documents/2020/03/environmental-cleaning-and-disinfectionprinciples-for-covid-19.pdf> Mar 18, 2020